

## If you adoted the following primer design:

A1 F1-R1

A2 F2-R1

A3 F3-R1

...

A4 F1-R2

A5 F2-R2

A6 F3-R2


...

A7 F1-R3

A8 F2-R3

A9 F3-R3

Your choice was the **combinatorial** strategy



**Figure 1:** With combinatorial dual indexing, every i5 and i7 is shared among other samples on the same plate with indexes repeating across rows and down columns; the combinations are unique though the individual indexes are not. In contrast to combinatorial dual indexing, unique dual indexing has completely unique indexes (eg, 96 unique i5s and 96 unique i7s per 96-well plate).

https://support.illumina.com/bulletins/2018/08/understanding-unique-dual-indexes--udi--and-associated-library-p.html

You can choose if you'd like toi merge the paired ends or to keep them separated, as needed for running DADA, for example, or for any other reason, but the combinatorial strategy is about tre primer design, not "merged or not". Right?


When you choose to merge your Illumina pairs before demultiplexing (for any reason), this is what happens in the merge step:



Mapping file: barcodes.txt

A1 CCTAAACTACGG CCATCACATAGG → BOTH FORWARD DIRECTION

Both pairs are in **forward** direction in the **original** fastq files, as you see above, but when they merge, the **R2 sequences become reverse** (see below, red). Because of this, when demultiplexing with linked mode, the vast majority of the merged reads are **TagF** – **reverse TagR** or **reverse TagF** – **TagR**:



Because of this I **recommend** to demultiplex **merged** pairs in **Linked adapter** mode, when both tags are **required** (and optionally **anchored**, in order to **minimize wrong** combinations during library preparation).